



FIRST RESPONSE[®]

HCV CARD TEST (Self Test) (Whole Blood)

For the detection of antibodies to Hepatitis C virus in human whole blood

One-step Hepatitis C Card Test for Whole blood

FEATURES


- ▶ Clear background.
- ▶ Self contained test procedure.
- ▶ Specifically engineered for whole blood specimen
- ▶ Storage temperature 4°C - 30°C.
- ▶ Detects all subtypes of HCV.
- ▶ Long Shelf life up to 24 months.

BENEFITS


- ▶ Unmistakable results at 15-20 minutes.
- ▶ Simple and easy to use, does not require advanced medical training.
- ▶ Cost effective, can be used effectively in resource-limited conditions.
- ▶ Accurate and reliable.
- ▶ Ideal for untrained/lay user.



RESULTS INTERPRETATION



TEST PROCEDURE



KIT CONTENTS

- Test device pouch
- Sterile safety lancet
- Alcohol swab
- Test Device
- Specimen transfer device
- Desiccant
- Dry swab
- Single use assay buffer vial
- Disposable biohazard bag
- Instructions for use
- Bandage

ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI03FRCST	10Self Tests / Kit



ISO 13485 & EN ISO 13485
Certified Company

Part No. : PI03-BRO-ST-001, Rev. AB, Date: 2021-11-11



Manufactured by

Premier Medical Corporation Private Limited

Manufacturing Site and Registered Office : A1-302, GIDC, Sarigam 396155. Dist. Valsad, Gujarat, INDIA

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,
• www.premiermedcorp.com